


Manifolds

529 SERIES

Protocol Purge Station

The 529 Series Protocol Purge Station combines all the safety and convenience of a standard Protocol Station with a feature to purge the hose and inlet side of the regulator. The Protocol Purge Station allows isolation of the regulator inlet and is equipped with a purge valve at the 12 o'clock position on the mounting block with a 1/4" stainless steel tube fitting to safely pipe away the process gas during purging, if required. This station comes complete with mounting block and bracket, two 533 Series diaphragm valves and 3-foot (900 mm) all stainless steel hoses with armor casing. Available in brass, chrome-plated brass, or 316L stainless steel.



Features

Regulator Isolation Valve

Isolates regulator inlet during purging process

Purge Valve with 1/4" Compression Tube Fitting

Allows complete and safe purging of the inlet

Integral Check Valve at Inlet of Hose

Minimizes purging during cylinder change

Bracket Mounts

Attaches conveniently to any surface

Connections

CGA, DIN 477, BS 341 and others available

Materials

Protocol Tee

Brass, chrome-plated brass or 316L stainless steel barstock

Flexible Hoses

316L stainless steel inner core Monel® inner core for oxygen

Inlet Connection

316L stainless steel or brass barstock

Check Valve "O" Ring

Viton® (or neoprene for ammonia)

Bracket

304 stainless steel

Specifications

Maximum Inlet Pressure

3500 PSIG (240 BAR)

Temperature Range

-40°F to 140°F (-40°C to 60°C)

Weight (Single Stage)

3.6 lbs. (1.4 kg)


Weight (Dual Stage)

4.0 lbs. (1.6 kg)

Manifolds



Installation Dimensions



529 0174 Dimensions Shown

Ordering Information		
Order Number	Description	Max Inlet
Add letter "Q" after any regulator order number.	For example, to order a 422 1331-580 with a protocol purge station, the order number would be 422 1331-580Q	
To order separately:		
529 0171-CON	Protocol purge station for single stage brass regulators	3500 PSIG (240 BAR)
529 0172-CON	Protocol purge station for single stage chrome-plated brass regulator	3500 PSIG (240 BAR)
529 0173-CON	Protocol purge station for single stage 316L stainless steel regulators	3500 PSIG (240 BAR)
529 0174-CON	Protocol purge station for dual stage brass regulators	3500 PSIG (240 BAR)
529 0175-CON	Protocol purge station for dual stage chrome-plated brass regulators	3500 PSIG (240 BAR)
529 0176-CON	Protocol purge station for dual stage 316L stainless steel regulators	3500 PSIG (240 BAR)