

405 SERIES REGULATORS



The 405 Series single stage regulators are intended for secondary pressure control of non-corrosive, ultra-high purity, or liquefied gases in high purity gas distribution systems. For applications where gas is piped up to instruments, the 405 Series point-of-use regulator allows the user to tightly monitor and adjust delivery pressure as needed.

Typical Applications

- Bulk Gas Distribution Systems
- High Purity Carrier Gases
- Control of Cryogenic Gases
- Gas and Liquid Chromatography
- Zero, Span, and Calibration Gases
- High Purity Chamber Pressurization
- Liquefied Hydrocarbon Gas Control



405 3031-01-000 shown

Features

Metal-to-Metal Diaphragm Seal prevents gas contamination

Brass Barstock Body provides smooth, durable surface finish

Cleanroom Assembly eliminates contamination at installation

Pipe Away Relief Valve vents exhaust gases safely

CAPSULE® Seat Design features integral 360° filtration in both stages

Pressure Ranges from 0-15 to 0-500 PSIG (0-1 to 0-34 BAR) accommodate a broad range of applications

Materials and Specifications

Maximum Inlet Pressure (bare body): 3000 PSIG (210 BAR)

Body: Brass barstock

Bonnet: Brass barstock

Gauge: 2 in (50 mm) diameter brass

Seat: PTFE

Filter: 10-micron sintered bronze

Temperature Range: -40°F to 140°F (-40°C to 60°C)

Diaphragm: 316L stainless steel

Cv: 0.1 Max outlet 50 PSIG (3.5 BAR) or below; 0.2 Max outlet above 50 PSIG (3.5 BAR) *See flow curves attached*

Helium Leak Integrity: 1×10^{-9} scc/sec

Internal Seals: PTFE

Ports (bare body): 1/4 in FNPT


Weight (405 3031-000): 2.71 lbs. (1.23 kg)

Conformances: Cleanliness meets or exceeds CGA G-4.1; PED 2014/68/EU; ANSI/ASME B40.1; CRN OH5216

405 SERIES REGULATORS



Installation Dimensions



Ordering Information

405	A	B	C	D	-CON	Options	
Series 405	Outlet Pressure	Outlet Gauge	Inlet Gauge	Outlet Assemblies	Assembly/Gauges	Inlet Connections	Installed Options
	1: 0-15 PSIG (0-1 BAR)	30"-0-30 PSIG/-1-0-2 BAR	0: None	0: 1/4" FNPT port	0: Bare body	000: 1/4" FNPT	D: Deep purge
	2: 0-50 PSIG (0-3.5 BAR)	30"-0-100 PSIG/-1-0-7 BAR		1: 1/4" MNPT	1: Cleanroom assembly (PSIG/kPa gauges)	TF2: 1/8" tube	M: Protocol station
	3: 0-100 PSIG (0-7 BAR)	30"-0-200 PSIG/-1-0-14 BAR		2: 1/4" tube fitting	2: Cleanroom assembly (BAR/PSIG gauges)	TF4: 1/4" tube	Q: Protocol purge station
	4: 0-250 PSIG (0-17 BAR)	0-400 PSIG/0-27 BAR		3: Diaphragm valve 1/4" tube fitting	6: Mirror image (PSIG/kPa gauges)	TF6: 3/8" tube	
	5: 0-500 PSIG (0-34 BAR)**	0-1000 PSIG/0-70 BAR		4: Diaphragm valve 1/4" MNPT	7: Mirror image (BAR/PSIG gauges)	M06: 6mm tube	
	7: 0-150 PSIG (0-10 BAR)	30"-0-200 PSIG/-1-0-14 BAR		5: Needle valve 1/4" MNPT			
	<i>*Standard assembly does not include relief valve.</i>			6: 1/8" tube fitting			
				7: 3/8" tube fitting			
				8: Diaphragm valve 1/8" fitting			
				9: Diaphragm valve 1/4" FNPT			
				A: 3/8" BSP RH fitting			
				B: Diaphragm valve 3/8" tube fitting			
				C: 3/8" BSP LH fitting			
				D: 6mm brass hose barb (Not available if A=4 or 5)			
				G: 1/8" stainless steel tube fitting			
				H: 1/4" stainless steel tube fitting			
				M: 6mm tube fitting			
				S: Diaphragm valve 6mm tube fitting			

*Standard assembly does not include relief valve.

Related Options

Part No.	Description
550 0001	Captured vent kit
550 0002	Panel mount kit
476 0002	Helium Leak certification

405 SERIES REGULATORS



Flow Curves for 405 Series

